

Safety Data Sheet

According to EC Regulation No. 453/2010/EC

Version No.: 7 Revision Date: 30 November 2023

Print Date: 30 Nov. 23

Section 1. Identification

.1 Product Identifier 00156208

Product Name CHLORIDE ANALYSER BUFFER

CAS Number Mixture
REACH Registration No Not applicable

1.2 Relevant identified uses of the substance or mixture & uses advised against

Uses of Material Chemical for industrial and laboratory use. Not suitable for domestic use.

1.3 Supplier Sherwood Scientific Ltd

1 The Paddocks Cherry Hinton Road

Cambridge CB1 8DH

UNITED KINGDOM

Phone +44(0)1223 243444

Email <u>hands@sherwood-scientific.com</u>

L.4 Emergency Telephone (08:00 – 16:30) 01223 243444

(24hr) 11

(Have this document to hand)

Section 2. Hazards Identification

2.1 Classification of the substance or mixture

Classification according to regulation 1272/2008/EC

Skin corrosion/irritation, category 2 (Skin Irrit. 2). Serious eye damage/irritation, category 2 (Eye Irrit. 2).

2.2 Label elements

Labelling according to regulation 1272/2008/EC

Signal word Warning

Hazard Pictograms



Hazard Statements

H315, Causes skin irritation., H319, Causes serious eye irritation.

Precautionary Statements

Wash thoroughly after handling. Wear protective gloves / protective clothing / eye protection. IF ON SKIN: Wash with soap and water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and

wash before reuse.

Section 3. Composition

3.2 Mixtures

Component	CAS No.	EEC No.	Conc w/w	CLP Classification (1272/2008/CE)
Acetic acid	64-19-7	200-580-7	10.0 %	Flam. Liq. 3,Skin Corr. 1A

Section 4. First Aid

4.1 Description of first aid measures

Irrigate thoroughly with plenty of water for at least 10 minutes, holding the eye open. OBTAIN MEDICAL Eyes

ATTENTION URGENTLY.

Skin Wash off skin thoroughly with water. Remove contaminated clothing immediately and wash before re-use.

If discomfort persists OBTAIN MEDICAL ATTENTION.

Inhalation Remove from exposure. Keep warm and at rest. If there is difficulty in breathing give oxygen if available. If

discomfort persists OBTAIN MEDICAL ATTENTION.

Ingestion If conscious give plenty of water to drink. Do not induce vomiting. If there is difficulty in breathing give oxygen

if available. If breathing stops or shows signs of failing, apply artificial resuscitation. If unconscious place in the recovery position. OBTAIN MEDICAL ATTENTION URGENTLY.

4.2 Most important symptoms and effects, both acute & delayed.

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed.

No further relevant information available.

Section 5. Fire Fighting

5.1 Extinguishing media

Extinguishing Media Water spray, alcohol resistant foam, dry powder or carbon dioxide. Use water spray to keep fire

exposed containers cool.

Unsuitable Media Do not use water jet.

5.2 Special hazards arising from the substance or mixture

Hazards Presents no specific fire danger.

5.3 Advice for firefighters

Advice for firefighters Evacuate area immediately. Keep up wind. Avoid exposure to toxic vapours and fumes. Fire-fighters should wear

protective clothing and breathing apparatus.

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal Protection Avoid breathing vapour. Use approved personal protective equipment. Evacuate area immediately. Do not

allow general use of area until it is safe to do so.

6.2 Environmental precautions

Environmental Keep material out of sewers, storm drains, surface waters and soil. Notify the Environmental Agency and local

Environmental Health Officer if major spillage occurs.

6.3 Methods and material for containment and cleaning up

Major Spillage Contain and absorb on inert material. Transfer absorbent to salvage container for removal. Wash area down

with copious amounts of water.

Minor Spillage Wash area down with copious amounts of water.

6.4 Reference to other sections

See section 8.2 for information on protective equipment and section 13 for information on disposal.

Section 7. Storage & Handling

7.1 Precautions for safe handling

Avoid contact with skin and eyes. Do not breath vapours. Do not allow to contaminate clothing.

Ensure Local Exhaust Ventilation maintains vapour concentrations below the recommended

limits.

7.2 Conditions for safe storage, including any incompatibilities

Well ventilated, cool, dry storage.

Storage Class: 12

7.3 Specific end use(s)

See section 1.2.

Section 8. Workplace Exposure & Personal Protection

8.1 Control parameters

Component	CAS No	Concentration	Workplace Exposure Limits				
			Long Term (8hr TWA)		Short Term (15min period)		
Acetic acid	64-19-7	10%	10.0 ppm	15.0 mg/m-3	25.0 ppm	37.0 mg/m-3	

8.2 Exposure controls

Respiratory Protection Use L.E.V. or natural ventilation to maintain vapour concentrations below exposure limits. If not, use a well

maintained chemical cartridge organic vapour respirator, or use self-contained breathing apparatus.

Hand Protection Use nitrile gloves or PVC gauntlets.

Eye Protection Use chemical full face shield.

Skin Protection Wear PVC oversuit.

Special Hazards No special precautions required.

Section 9. Physical & Chemical Properties

9.1 Information on basic physical and chemical properties

Appearance Clear colourless liquid.

Odour Sharp vinegary odour and burning taste. 5 @ 20°C

Boiling Point Not available
Melting Point Not applicable
Flash Point Not applicable
Upper Flammable Limit Not applicable
Lower Flammable Limit Not applicable
Auto Ignition Not applicable

Explosive Properties No.

Oxidising Properties No.

Vapour Pressure Not applicable Relative Density 1.1160

9.2 Other information

No data available.

Section 10. Stability & Reactivity

10.1 Reactivity No data available.

10.2 Chemical Stability Stable under normal conditions

10.3 Possibility of hazardous

reactions

No data available.

10.4 Conditions to Avoid No specific conditions.

10.5 Incompatible Materials Hydrogen peroxide, chromium trioxide and potassium permanganate. Potassium t-butoxide. Alkalis.

10.6 Hazardous Decomposition

Products

None unusual. Burning will produce smoke, carbon monoxide and/or carbon dioxide

Section 11. Toxicological Information

11.1 Information on toxicological effects

Eyes Both the vapour and liquid will, be irritating to the eyes. Skin Both the vapour and liquid will, be irritating to the skin.

LD50 Skin Not available

Causes severe corrosion of the mouth, throat and gastro-intestinal Ingestion

tract. LD50 Oral

Inhalation High concentrations of vapour will produce severe irritation of the eyes, nose, throat and respiratory tract.

LD50 Inhalation Not available **TCLo** Not available

Carcinogenicity Not considered to be a carcinogen. Mutagenicity Not considered to be a mutagen.

Reproductive Effects None identified.

Other Information The irritant effect provides warning that control of exposure is needed. 10ppm is the threshold for irritation

with severe irritation occurring above 25ppm.

Section 12. Ecological

12.1 Toxicity Readily bio-degraded in the environment. Slightly toxic to aquatic species but is unlikely to bioaccumulate.

LC50 Algal Not available LC50 Crustacea Not available LC50 Fish Not available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential No data available. 12.4 Mobility in soil No data available.

12.5 Results of PBT & vPvB

Assessment

Assessment not required.

12.6 Other adverse effects None known at present:

Section 13. Disposal Considerations

13.1 Waste treatment methods

Disposal Methods Dilute in a large excess of water and carefully neutralise with soda ash, then wash to drain with copious amounts

of water.

Contaminated Packaging Wash out containers with water. Use a licensed waste disposer.

Section 14. Transport Information

14.1 UN Number Non-restricted14.2 Proper Shipping Name Non-restricted

14.3 Transport classes

UN classification None Subsidiary Risk None Transport category None

ADR Hazard ID Non-restricted
Tunnel Restriction Code Non-restricted

14.4 Packing Group None

14.5 Environment hazards See section 12.

14.6 Special precautions for No special precautions required.

user

14.7 Transport in bulk Not transported in bulk.

Section 15. Regulatory Information

15.1 Safety, health and environment regulations specific for substance/mixture.

Classification, Labeling & Packaging of Substances & Mixtures Regulations (1272/2008/CE)

Signal word Warning

Hazard Pictograms



Hazard Statements H315, H319

Causes skin irritation. Causes serious eye irritation.

German Water Hazard Class: WGK 1

Precautionary Statements P264, P280, P302+P352, P332+P313, P362

Wash thoroughly after handling. Wear protective gloves / protective clothing / eye protection. IF ON SKIN: Wash with soap and water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse.

Chemicals (Hazard Information & Packaging) Regulations 2009 (67/548/EEC, 99/45/EC)

Classification Irritant (Xi)



Risk & Safety Phrases R36/38, S24/25, S26

Irritating to eyes and skin. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

15.2 Chemical safety assessment

Assessment not required.

Section 16. Other Information

The information contained in this document only covers the hazards presented by this material, it DOES NOT constitute a workplace risk assessment. See sections 11 for toxicological information and section 12 for ecological information.